

## SA-300VMA

Vertical Autoclave With Vapor  
50 Liter



- CLASS-B Autoclave
- Digital display
- LED display for program (function)
- Sterilization temperature: 121°C / 135°C
- Programs for selection 121°C / 135°C, wrapped / unwrapped, dry / undry, Liquid, PRION program, Bowie-deck, Pump test
- Indication: Door closed, Low water, Over heat, Over pressure...
- Auto dry function, dry time can be selected
- External dimensions (dxwxh): 480mm x 650mm x 1084mm
- Chamber size: 300mm x 700mm
- Chamber diameter: 300mm
- Analog vacuum gauge
- Auto-Fill water system
- 2 Baskets

Display of SA-300VMA



Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11  
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04  
e-mail: info@uniequip.de · Internet: www.uniequip.com

UNI-EQUIP

## Horizontal Table Top Autoclave



- Pre-Vacuum & Dry-Vacuum
- CLASS B – TÜV and CE certified, ISO 9001
- 16 or 24 liters chamber capacity
- Sterilization temperature 121°C / 135°C
- Intelligent design for your professional choice
- Installed 24 programs and special programs for PRION sterilization, LIQUID sterilization and test programs, full automatic processing
- Auto-Fill water system
- Door lock indication
- Digital temperature display
- Analog vacuum gauge

UNI-EQUIP

Model	SA-260MA	SA-230MA
Auto dry time can be selected. 0, 15, 30, 35, 40 mins	Yes	Yes
Ext dimensions (dxwxh)	655mm x 533mm x 442mm	594mm x 522mm x 410mm
Water capacity of water tank	4200cc	4200cc
Chamber capacity	24 liters	16 liters
Net weight	52 kg	49 kg
Power supply	230V 50/60Hz	230V 50/60Hz
Power consumption	2400W / 11Amp.	1900W / 8Amp.
<b>SAFETY DEVICE</b>		
Pressure safety valve & protection switch	Yes	Yes
Over temperature protection	Yes	Yes
Water tank water level sensor	Yes	Yes
Pressure door auto-lock	Yes	Yes
Standard accessories	Sterilization Box, Tray 1 pc.	Sterilization Box, Tray 1 pc.
<b>OPTIONAL ACCESSORIES</b>	Sterilization Plate set	Sterilization Plate set

Dimensions and specifications are subject to change without notice for further improvements on design.  
Dimensions may vary dependent on different accessories be used.

UNI-EQUIP

D-82152 Martinsried / München - Fraunhoferstrasse 11  
Telefon: +49-89-8 57 52 00 · Telefax: +49-89-8 56 13 04  
e-mail: info@uniequip.de · Internet: www.uniequip.com