

Urine Adulteration Tests

Drug Adulteration Strip, *colorimetric test*

DOA Saliva Tests

CE

Alcohol Saliva Strip

DOA Saliva Multi-6 Test

(AMP, COC, MET, MOP, PCP, THC)

DOA Saliva Multi-6/1 Test

(AMP, COC, MET, MOP, MTD, THC)

both also available as High Sensitivity Touch Test

Pregnancy / Fertility

CE

hCG Midstream *(urine)**

hCG Midstream for Home Use *(urine)**

hCG Dipstick *(urine)*

hCG Combo Dipstick *(urine/serum)*

hCG Dipstick for Home Use *(urine)*

hCG Cassette *(urine)*

hCG Combo Cassette *(urine/serum)*

LH Dipstick *(urine)*

LH Dipstick for Home Use *(urine)*

LH Cassette *(urine)*

LH Cassette for Home Use *(urine)*

** upon request*

Tumor Markers

AFP Cassette *(serum/plasma/whole blood)* **CE**

CEA Cassette *(serum/plasma/whole blood)*

FOB Cassette *(human fecal occult blood)* **CE**

FOB-TRF Combi Cassette *(stool samples)* **CE** **new**

PSA Cassette *(serum/plasma)*

Cardiac Markers

CE

Cardiac Combo Cassette *(serum/plasma/whole blood)*

Troponin I Cassette *(serum/plasma/whole blood)*

Others


CE

CRP Dipstick *(whole blood)*

D-Dimer Cassette *(plasma/whole blood)*

Microalbumin Dipstick *(urine)*

>>> **DIAQUICK RAPID TESTS**



Aids, Hepatitis

HAV IgM Cassette (serum/plasma/whole blood) **new**

HCV Ab Cassette (serum/plasma)

HCV Ab Cassette (serum/plasma/whole blood)

HBsAb Cassette (serum/plasma)

HBsAb Dipstick (serum/plasma)

HIV 1&2 Ab Cassette (serum/plasma/whole blood) **CE**

HIV Triline Cassette (serum/plasma/whole blood)

HIV Ab Cassette (serum/plasma)

HIV Ab Cassette (serum/plasma/whole blood)

HBsAg Cassette (serum/plasma)

HBsAg Dipstick (serum/plasma)

HBsAg Cassette (serum/plasma/whole blood)

Other Infectious Diseases **CE**

Adenovirus Cassette (stool samples)

Chikungunya IgM Cassette (serum/plasma/whole blood)*

Chlamydia Cassette (endocervical/urethral/urine samples)

Dengue IgG/IgM Ab Cassette (serum/plasma/whole blood)

H.pylori Cassette (serum/plasma/whole blood)

H.pylori Stool Cassette (stool samples)

Influenza Ag Dipstick (nasal/throat swab samples)

Influenza Ag A/B/H1N1 Dipstick (nasal/throat swab samples)*

Malaria P.f. Cassette (whole blood)

Malaria P.f./Pan Cassette (whole blood)

Mononucleosis Cassette (serum/plasma/whole blood)

Mycoplasma IgM Cassette (serum/whole blood) **new**

Rotavirus Cassette (stool samples)

Rota/Adeno Combi Cassette (stool samples)

S.typhi IgG/IgM Ab Cassette (serum/plasma/whole blood)**

S.typhi Ag Cassette (stool/serum/plasma/whole blood)**

S.typhi/paratyphi Ag Cassette (stool/serum/plasma/whole blood)**

Strep.A Cassette (throat swab samples)

Strep.A Dipstick (throat swab samples)

Syphilis Ab Cassette (serum/plasma/whole blood)

Syphilis Ab Dipstick (serum/plasma/whole blood)

Tuberculosis IgG/IgM Cassette (serum/plasma)*, **

* upon request ** no **CE**