

EPi® Drawer Racks

Width 140 mm. Stainless steel.


With Folding Handle Open Design



Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

With Folding Handle Closed Design




Please add to cat. no.:

- L** Base of drawer in middle open.
Facilitates removal and insertion.
- V** Base of drawer closed

EPi® 45 + 50 Drawer Racks. With Folding Handle

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Front closed.

With Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33NAM	■	3	3	434 x 184
SE50-34NAM	■	3	4	434 x 245
SE50-35NAM	■	3	5	434 x 306
SE50-36NAM	■	3	6	434 x 367
SE50-43NAM	■	4	3	567 x 184
SE50-44NAM	■	4	4	567 x 245
SE50-45NAM	■	4	5	567 x 306
SE50-46NAM	■	4	6	567 x 367
SE50-53NAM	■	5	3	700 x 184
SE50-54NAM	■	5	4	700 x 245
SE50-55NAM	■	5	5	700 x 306
SE50-56NAM	■	5	6	700 x 367




Without Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33NAO	■	3	3	416 x 173
SE50-34NAO	■	3	4	416 x 231
SE50-35NAO	■	3	5	416 x 288
SE50-36NAO	■	3	6	416 x 346
SE50-43NAO	■	4	3	549 x 173
SE50-44NAO	■	4	4	549 x 231
SE50-45NAO	■	4	5	549 x 288
SE50-46NAO	■	4	6	549 x 346
SE50-53NAO	■	5	3	682 x 173
SE50-54NAO	■	5	4	682 x 231
SE50-55NAO	■	5	5	682 x 288
SE50-56NAO	■	5	6	682 x 346




EPi® 45 + 50 Drawer Racks. With Folding Handle

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Front closed.

With Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33GAM	■	3	3	444 x 184
SE50-34GAM	■	3	4	444 x 245
SE50-35GAM	■	3	5	444 x 306
SE50-36GAM	■	3	6	444 x 367
SE50-43GAM	■	4	3	577 x 184
SE50-44GAM	■	4	4	577 x 245
SE50-45GAM	■	4	5	577 x 306
SE50-46GAM	■	4	6	577 x 367
SE50-53GAM	■	5	3	710 x 184
SE50-54GAM	■	5	4	710 x 245
SE50-55GAM	■	5	5	710 x 306
SE50-56GAM	■	5	6	710 x 367



Without Safety Stop				
Cat. No.	L	Number	DxH	
	V	Depth	Height	(mm)
SE50-33GAO	■	3	3	426 x 173
SE50-34GAO	■	3	4	426 x 231
SE50-35GAO	■	3	5	426 x 288
SE50-36GAO	■	3	6	426 x 346
SE50-43GAO	■	4	3	559 x 173
SE50-44GAO	■	4	4	559 x 231
SE50-45GAO	■	4	5	559 x 288
SE50-46GAO	■	4	6	559 x 346
SE50-53GAO	■	5	3	692 x 173
SE50-54GAO	■	5	4	692 x 231
SE50-55GAO	■	5	5	692 x 288
SE50-56GAO	■	5	6	692 x 346




With Hand Grip Open Design



With Hand Grip Closed Design



Without Handle (Tray) Without Safety Stop



EPi® 45 + 50 Drawer Racks. With Hand Grip

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

With Safety Stop			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33NBMV	3	3	434 x 184
SE50-34NBMV	3	4	434 x 245
SE50-35NBMV	3	5	434 x 306
SE50-36NBMV	3	6	434 x 367
SE50-43NBMV	4	3	567 x 184
SE50-44NBMV	4	4	567 x 245
SE50-45NBMV	4	5	567 x 306
SE50-46NBMV	4	6	567 x 367
SE50-53NBMV	5	3	700 x 184
SE50-54NBMV	5	4	700 x 245
SE50-55NBMV	5	5	700 x 306
SE50-56NBMV	5	6	700 x 367

Without Safety Stop			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33NB0V	3	3	434 x 173
SE50-34NB0V	3	4	434 x 231
SE50-35NB0V	3	5	434 x 288
SE50-36NB0V	3	6	434 x 346
SE50-43NB0V	4	3	567 x 173
SE50-44NB0V	4	4	567 x 231
SE50-45NB0V	4	5	567 x 288
SE50-46NB0V	4	6	567 x 346
SE50-53NB0V	5	3	700 x 173
SE50-54NB0V	5	4	700 x 231
SE50-55NB0V	5	5	700 x 288
SE50-56NB0V	5	6	700 x 346

EPi® 45 + 50 Drawer Racks. With Hand Grip

Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

With Safety Stop			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33GBMV	3	3	444 x 184
SE50-34GBMV	3	4	444 x 245
SE50-35GBMV	3	5	444 x 306
SE50-36GBMV	3	6	444 x 367
SE50-43GBMV	4	3	577 x 184
SE50-44GBMV	4	4	577 x 245
SE50-45GBMV	4	5	577 x 306
SE50-46GBMV	4	6	577 x 367
SE50-53GBMV	5	3	710 x 184
SE50-54GBMV	5	4	710 x 245
SE50-55GBMV	5	5	710 x 306
SE50-56GBMV	5	6	710 x 367

Without Safety Stop			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33GB0V	3	3	444 x 173
SE50-34GB0V	3	4	444 x 231
SE50-35GB0V	3	5	444 x 288
SE50-36GB0V	3	6	444 x 346
SE50-43GB0V	4	3	577 x 173
SE50-44GB0V	4	4	577 x 231
SE50-45GB0V	4	5	577 x 288
SE50-46GB0V	4	6	577 x 346
SE50-53GB0V	5	3	710 x 173
SE50-54GB0V	5	4	710 x 231
SE50-55GB0V	5	5	710 x 288
SE50-56GB0V	5	6	710 x 346

EPi® 45 + 50 Drawer Racks. Without Handle


Useful for EPi® 45 + 50 and all boxes with external dimensions of up to 133x133x53 mm (WxDxH). Base of drawers closed.

Open Design			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33NTOV	3	3	416 x 173
SE50-34NTOV	3	4	416 x 231
SE50-35NTOV	3	5	416 x 288
SE50-36NTOV	3	6	416 x 346
SE50-43NTOV	4	3	549 x 173
SE50-44NTOV	4	4	549 x 231
SE50-45NTOV	4	5	549 x 288
SE50-46NTOV	4	6	549 x 346
SE50-53NTOV	5	3	682 x 173
SE50-54NTOV	5	4	682 x 231
SE50-55NTOV	5	5	682 x 288
SE50-56NTOV	5	6	682 x 346

Closed Design			
Cat. No.	Number		DxH (mm)
	Depth	Height	
SE50-33GTOV	3	3	426 x 173
SE50-34GTOV	3	4	426 x 231
SE50-35GTOV	3	5	426 x 288
SE50-36GTOV	3	6	426 x 346
SE50-43GTOV	4	3	559 x 173
SE50-44GTOV	4	4	559 x 231
SE50-45GTOV	4	5	559 x 288
SE50-46GTOV	4	6	559 x 346
SE50-53GTOV	5	3	692 x 173
SE50-54GTOV	5	4	692 x 231
SE50-55GTOV	5	5	692 x 288
SE50-56GTOV	5	6	692 x 346

EPPI® 32 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. SE32-38NCMV
For 24 x EPPI® 32 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x325 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 37 Drawer Racks

With Grip Rail Open Design




Example: Cat. No. SE37-46NCML
For 24 x EPPI® 37 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x274 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 32 Drawer Racks


Without Handle (Tray) Open Design



Example: Cat. No. SE32-28NTOV
For 16 x EPPI® 32 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x283x297 mm


EPPI® 37 Drawer Racks

Without Handle (Tray) Open Design



Example: Cat. No. SE37-35NTOV
For 15 x EPPI® 37 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x211 mm

With Grip Rail




Without Handle (Tray)




EPPi® 75 Drawer Racks

With Folding Handle Open Design¹




EPPi® 75 Drawer Racks

With Hand Grip Open Design¹




EPPi® 75 Drawer Racks

Without Handle (Tray) Open Design¹




EPPi® 80 Drawer Racks

With Folding Handle Open Design¹




EPPi® 80 Drawer Racks

With Hand Grip Open Design¹



EPPi® 80 Drawer Racks


Without Handle (Tray) Open Design¹



¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33).
Please enquire if needed.

EPPI® 96 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE96-43NAOL
For 12 x EPPI® 96 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x549x305 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 96 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE96-33NBMV
For 9 x EPPI® 96 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x316 mm

EPPI® 96 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SE96-23NTOV
For 6 x EPPI® 96 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x283x305 mm

EPPI® 102 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE102-32NAML
For 6 x EPPI® 102 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x434x223 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 102 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE102-52NBVO
For 10 x EPPI® 102 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x700x216 mm

EPPI® 102 Drawer Racks.

Without Handle (Tray) Open Design¹




Example: Cat. No. SE102-42NTOV
For 8 x EPPI® 102 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x549x216 mm



- ¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33). Please enquire if needed.
- ² All drawer racks can be produced with additional drawers for shallower boxes (e.g. EPPI® 50) to exploit the interior space of uprights to the maximum.

EPPI® 105 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE105-43NAML
For 12 x EPPI® 105 Boxes.
With safety stop. Base of drawers open.
Dimensions: WxDxH = 140x567x340 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 105 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE105-32NBMV
For 6 x EPPI® 105 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x434x227 mm

EPPI® 105 Drawer Racks


Without Handle (Tray) Open Design¹



Example: Cat. No. SE105-53NTOV
For 15 x EPPI® 105 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x682x329 mm

EPPI® 128 Drawer Racks

With Folding Handle Open Design¹




Example: Cat. No. SE128-22NAOL
For 4 x EPPI® 128 Boxes.
Without safety stop. Base of drawers open.
Dimensions: WxDxH = 140x283x267 mm

L Base of drawer in middle open. Facilitates removal and insertion.
V Base of drawer closed.

See illustrations on page 32

EPPI® 128 Drawer Racks


With Hand Grip Open Design¹



Example: Cat. No. SE128-42NBMV
For 8 x EPPI® 128 Boxes.
With safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x567x274 mm

EPPI® 128 Drawer Racks.

Without Handle (Tray) Open Design¹



Example: Cat. No. SE128-32NTOV
For 6 x EPPI® 128 Boxes.
Without safety stop. Base of drawers closed.
Dimensions: WxDxH = 140x416x267 mm



¹ All drawer racks can also be supplied with a closed design. (See illustrations at page 32+33). Please enquire if needed.

² All drawer racks can be produced with additional drawers for shallower boxes (e.g. EPPI® 50) to exploit the interior space of uprights to the maximum.

EPPI® Upright Bins


Width 139 and 278 mm. Stainless steel.



EPPI® Upright Bins

Single Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK169-121 BN	1	168	139x278x169
SEK169-131 BN	1	168	139x414x169
SEK169-141 BN	1	168	139x549x169
SEK169-151 BN	1	168	139x685x169
SEK224-121 BN	1	223	139x278x224
SEK224-131 BN	1	223	139x414x224
SEK224-141 BN	1	223	139x549x224
SEK224-151 BN	1	223	139x685x224
SEK250-121 BN	1	249	139x278x250
SEK250-131 BN	1	249	139x414x250
SEK250-141 BN	1	249	139x549x250
SEK250-151 BN	1	249	139x685x250


Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK279-121 BN	1	278	139x278x279
SEK279-131 BN	1	278	139x414x279
SEK279-141 BN	1	278	139x549x279
SEK279-151 BN	1	278	139x685x279
SEK334-121 BN	1	333	139x278x334
SEK334-131 BN	1	333	139x414x334
SEK334-141 BN	1	333	139x549x334
SEK334-151 BN	1	333	139x685x334
SEK389-121 BN	1	388	139x278x389
SEK389-131 BN	1	388	139x414x389
SEK389-141 BN	1	388	139x549x389
SEK389-151 BN	1	388	139x685x389

EPPI® Upright Bins

Double Width

With same external dimensions as horizontal racks. Open at top. With hand grip. Ideal for combining with horizontal racks and drawer racks.



Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK169-221 BN	1	168	278x278x169
SEK169-231 BN	1	168	278x414x169
SEK169-241 BN	1	168	278x549x169
SEK169-251 BN	1	168	278x685x169
SEK224-221 BN	1	223	278x278x224
SEK224-231 BN	1	223	278x414x224
SEK224-241 BN	1	223	278x549x224
SEK224-251 BN	1	223	278x685x224
SEK250-221 BN	1	249	278x278x250
SEK250-231 BN	1	249	278x414x250
SEK250-241 BN	1	249	278x549x250
SEK250-251 BN	1	249	278x685x250

Standard Heights

Cat. No.	Number of Compartments	Inner Height mm	Ext. Dimensions WxDxH (mm)
SEK279-221 BN	1	278	278x278x279
SEK279-231 BN	1	278	278x414x279
SEK279-241 BN	1	278	278x549x279
SEK279-251 BN	1	278	278x685x279
SEK334-221 BN	1	333	278x278x334
SEK334-231 BN	1	333	278x414x334
SEK334-241 BN	1	333	278x549x334
SEK334-251 BN	1	333	278x685x334
SEK389-221 BN	1	388	278x278x389
SEK389-231 BN	1	388	278x414x389
SEK389-241 BN	1	388	278x549x389
SEK389-251 BN	1	388	278x685x389