



## ALPHA Box 136 x 136 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 136 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


| Colour |    |                         | Heights |       |       |       |        |        |        |
|--------|----|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| 1 or 2 |    |                         | 25 mm   | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| white  | WH | 1 standard<br>2 special | 025     | 032   | 050   | 075   | 100    | 110    | 130    |
| blue   | BU | 1 standard<br>2 special | 025     | 032   | 050   | 075   | 100    | 110    | 130    |
| yellow | YE | 1 standard<br>2 special | 025     | 032   | 050   | 075   | 100    | 110    | 130    |
| green  | GR | 1 standard<br>2 special | 025     | 032   | 050   | 075   | 100    | 110    | 130    |
| red    | RE | 1 standard<br>2 special | 025     | 032   | 050   | 075   | 100    | 110    | 130    |

How to determine  
the item no.:



Example 50 mm height,  
white, standard:

**KA 1 050 WH 1**



## Grid Dividers for ALPHA Box

|                 |       |               | Heights     |             |             |             |             |             |
|-----------------|-------|---------------|-------------|-------------|-------------|-------------|-------------|-------------|
|                 |       |               | 15 mm       | 25 mm       | 30 mm       | 40 mm       | 65 mm       | 80 mm       |
| Cells for vials |       |               | Cat. No.    | Cat. No.    | Cat. No.    | Cat. No.    | Cat. No.    | Cat. No.    |
| 16x16           | up to | ca. 7.5 mm Ø  | RA10151616N | RA10251616N | RA10301616N | RA10401616N | RA10651616N | RA10801616N |
| 15x15           | up to | ca. 8.0 mm Ø  | RA10151515N | RA10251515N | RA10301515N | RA10401515N | RA10651515N | RA10801515N |
| 14x14           | up to | ca. 8.7 mm Ø  | RA10151414N | RA10251414N | RA10301414N | RA10401414N | RA10651414N | RA10801414N |
| 13x13           | up to | ca. 9.4 mm Ø  | RA10151313N | RA10251313N | RA10301313N | RA10401313N | RA10651313N | RA10801313N |
| 12x12           | up to | ca. 10.2 mm Ø | RA10151212N | RA10251212N | RA10301212N | RA10401212N | RA10651212N | RA10801212N |
| 11x11           | up to | ca. 11.2 mm Ø | RA10151111N | RA10251111N | RA10301111N | RA10401111N | RA10651111N | RA10801111N |
| 10x10           | up to | ca. 12.5 mm Ø | RA10151010N | RA10251010N | RA10301010N | RA10401010N | RA10651010N | RA10801010N |
| 9x9             | up to | ca. 13.9 mm Ø | RA10150909N | RA10250909N | RA10300909N | RA10400909N | RA10650909N | RA10800909N |
| 8x8             | up to | ca. 15.7 mm Ø | RA10150808N | RA10250808N | RA10300808N | RA10400808N | RA10650808N | RA10800808N |
| 7x7             | up to | ca. 17.9 mm Ø | RA10150707N | RA10250707N | RA10300707N | RA10400707N | RA10650707N | RA10800707N |
| 6x6             | up to | ca. 21.0 mm Ø | RA10150606N | RA10250606N | RA10300606N | RA10400606N | RA10650606N | RA10800606N |
| 5x5             | up to | ca. 25.5 mm Ø | RA10150505N | RA10250505N | RA10300505N | RA10400505N | RA10650505N | RA10800505N |
| 4x4             | up to | ca. 32.0 mm Ø | RA10150404N | RA10250404N | RA10300404N | RA10400404N | RA10650404N | RA10800404N |
| 3x3             | up to | ca. 42.9 mm Ø | RA10150303N | RA10250303N | RA10300303N | RA10400303N | RA10650303N | RA10800303N |

When used with item number given in the overview (beginning with "RA1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x10 cells.

**KA1050WH1 + 0301010N = Cat. No. KA1050WH10301010N**