


## Kryo Box 100 / 10x10 • height 52 mm fix.

100-place cryogenic storage box made of extra strong polycarbonate.  
For 10x10 cryo vials 1.0/2.0 ml with female thread. Transparent lid with numeric ID code (1-100). 133x133x52 mm (WxDxH).

Cat. No.	Colour	Height
KPC50-C100BU	blue	52 mm
KPC50-C100GR	green	52 mm
KPC50-C100RE	red	52 mm
KPC50-C100YE	yellow	52 mm



## Kryo Box 81 / 9x9 • height 52 mm fix.

81-place cryogenic storage box made of extra strong polycarbonate.  
For 9x9 cryo vials 1.0/2.0 ml. Transparent lid with numeric ID code (1-81). 133x133x52 mm (WxDxH).


Cat. No.	Colour	Height
KPC50-C81BU	blue	52 mm
KPC50-C81GR	green	52 mm
KPC50-C81RE	red	52 mm
KPC50-C81YE	yellow	52 mm



## Kryo Box 25 / 5x5 • height 52 mm fix.

25-place cryogenic storage box made of extra strong polycarbonate.  
For 5x5 cryo vials 1.0/2.0 ml. Transparent lid with numeric ID code (1-25). 76x76x52 mm (WxDxH).

Cat. No.	Colour	Height
KPC50-C25BU	blue	52 mm
KPC50-C25GR	green	52 mm
KPC50-C25RE	red	52 mm
KPC50-C25YE	yellow	52 mm



## Kryo Box 81H / 9x9 • height 95 mm fix.

81-place cryogenic storage box made of extra strong polycarbonate.  
For 9x9 cryo vials 3.0 up to 5.0 ml. Transparent lid with numeric ID code (1-81). 133x133x95 mm (WxDxH).

Cat. No.	Colour	Height
KPC95-C81BU	blue	95 mm
KPC95-C81GR	green	95 mm
KPC95-C81RE	red	95 mm
KPC95-C81YE	yellow	95 mm

Matching Racks

Vertical Racks  
page 29

Horizontal Racks  
page 30-31

Drawer Racks  
page 32-33

Matching Racks

Vertical Racks  
page 29

Horizontal Racks  
page 30-31

Drawer Racks  
page 32-33

Matching Racks

Vertical Racks  
page 64

Horizontal Racks  
page 64

Drawer Racks  
Please request!

Matching Racks

Vertical Racks  
page 29

Horizontal Racks  
page 30-31

Drawer Racks  
page 36