

REMEMBER TO GET
A QUOTE ON
OPTIONS

MULTIPLE
RUNNING
TEMPERATURES

-85/-40°C

UPRIGHT
ULTRA LOW
TEMPERATURE
FREEZERS

ULUF 850

ULUF (Ultra Low Upright Freezer)

With its massive 826L the new ULUF 850 is still the most versatile Ultra low freezer in the market. From the possibility to disassemble for easy installation, to multiple temperature ranges at once, Arctiko set the standards in the Ultra low temperature freezing.



Racks and boxes.
Optional.




Multi flexible drawer
system for sample
storage. Optional.



Micro processor
controller with
digital display



Prepared for our
advanced CO2
backup system
(optional).



30 racks and 504
pieces 50mm plastic
boxes for the ULUF
850 (optional).

Model	ULUF 400	ULUF 450	ULUF 800	ULUF 850
Item code (DAI no.)	DAI 1407	DAI 1408	DAI 1409	DAI 1411
Dimensions exterior WxDxH (mm)	720x885x2089	720x885x2089	1492x885x2089	1492x885x2089
Net capacity (L)	413	413	826	826
Temperature range (°C)	-40/-20	-85/-40	-40/-20	-85/-40
Power supply (V)	230	230	230	230
Frequency (Hz)	50/60	50/60	50/60	50/60

We reserve the right to change specifications without notice. Subject to confirmation and availability.

Scan me and see more
about ULUF 850



www.arctiko.com